## DRY BATH -1No.

(Tender Ref. No.NIN/ST/12/Tender/2019/Dry bath/18)

## Features:

- LED display for current temperature and count down time
- Metal Blocks for easy replacement, cleaning and disinfecting
- Over temperature alarm protection
- Transparent cover to reduce the experiment error
- Self-Detection and alarm system

## **Specifications:**

- Cooling/heating setting range : -10°C to 100°C with very rapid cool down and heat up times
- Block dimensions: 100 x 110mm (Approx)
- Temperature range : -10°C to 100°C
- Stability: ± 0.1°C
- Display: LCD
- Heat up time: 25°C to 37°C ( 3-5mins)
  25°C to 100°C (15 25mins)
- Cool down time: 100°C to -10°C (20 -30mins)
  25°C to -10°C (15 25mins)
- Capacity: 20 x 0.5ml, 20x1.5ml, 20x2ml (micro tubes)
- Timer: 1min to 96 hours or non-stop
- Power supply: 230V AC, 50Hz